LICATIONS

Instant Sensitivity Gains Without the Need to Dry Down and Reconstitute Using Microelution Solid Phase Extraction (SPE)

Shahana Huq and Matt Brusius Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

Solid Phase Extraction (SPE) is an excellent cleanup solution for bioanalytical samples because it is selective, reproducible, and results in ultra-clean and concentrated samples. While the technique is effective, traditional formats such as 10 mg 96-well plates or tubes pose challenges for small or limited sample volumes and can result in dilute eluents if the sample is not dried down and reconstituted after the extraction procedure. This dry down step can require 30 or more minutes, adding a significant amount of time to the procedure. To overcome these challenges, the microelution 96-well SPE plate format was developed. The microelution SPE plates contain significantly less sorbent as compared to a traditional 10 mg 96-well SPE plate, allowing analysts to elute in volumes as low as 25 µL. These low elution volumes result in ultra-concentrated samples that do not need to be dried down. In this study, we compared the microelution 96-well SPE plate to a traditional 10 mg 96-well SPE plate to determine the amount of sensitivity gain that can be achieved by moving a 10 mg SPE method to the microelution format.

Experimental Conditions

Extraction Procedures

Timolol and procaine were extracted from serum using a traditional 10 mg 96-well SPE plate and a 2 mg microelution 96-well SPE plate, each of which was packed with the same polymeric strong cation-exchange chemistry, Strata®-X-C. The same volume of sample was loaded onto each SPE plate and the condition, equilibration, wash, and elution solvents were the same across both extractions to standardize the extraction protocols.

Table 1.

SPE Extraction of Timolol and Procaine from Serum

	Strata-X-C 96-Well SPE Plate, 10 mg/well Strata-X-C Microelution 96-Well SPE Plate, 2 mg/well			
Condition	500 µL Methanol	200 µL Methanol		
Equilibrate	500 µL Water	200 µL Water		
Load	$750\mu L$ diluted serum (375 μL serum diuted 1:1 with 4 % Phosphoric acid in water)	$750\mu\text{L}$ diluted serum (375 μL serum diuted 1:1 with 4 % Phosphoric acid in water)		
Wash 1	500 µL 2 % Formic acid in water	200 µL 2 % Formic acid in water		
Wash 2	500 µL Methanol	200 µL Methanol		
Elute	375 μL (3x 125 μL) 5 % Ammonium hydroxide in acetonitrile/methanol (60:40)	25 µL 5 % Ammonium hydroxide in acetonitrile/methanol (60:40)		
Inject	1µL	1 μL		



Matt Brusius Product Manager, Sample Preparation Matt Brusius is an avid ice hockey player. He likes skating



APPLICATIONS



HPLC Conditions

After cleanup by SPE, 1 μ L of each extraction was injected onto a Kinetex[®] 2.6 μ m Biphenyl core-shell HPLC column and the resulting concentrations were determined by LC/MS/MS.

Column: Dimensions: Part No.: Mobile Phase:	•••		
	0.00 5		
Gradient:	Time (min)	B (%)	
	0.00	5	
	0.50	5	
	3.00	95	
	3.50	95	
	3.51	5	
	5.50	5	
Flow Rate:	500 µL/min		
Temperature:	Ambient		
Detection:	MS/MS, API	5000™ (AB SCIEX)	

Results and Discussion

Both timolol and procaine were extracted from serum using two different SPE formats; a traditional 10 mg 96-well SPE plate and a 2 mg microelution 96-well SPE plate. The microelution format is designed to allow analysts to elute in small sample volumes which results in ultra-concentrated samples without the need to perform a dry down step. To standardize our comparison, the same polymeric strong cation-exchange sorbent, Strata®-X-C, was packed in each 96-well SPE plate and the same method was performed however the solvent volumes were reduced for the microelution method (**Table 1**). While both extraction methods were extremely effective and resulted in clean samples, the microelution format resulted in 15x more concentrated samples (**Figures 1** and **2**) when 1μ L of the eluent was injected onto the LC/MS/MS.

Figure 1. Timolol (1ng/mL) Extracted from Serum Resulted in 15x More Concentrated Samples using the microelution Format





Figure 2. Procaine (1ng/mL) Extracted from Serum Resulted in 15x More Concentrated Samples using the microelution Format



Conclusion

While both the 10 mg 96-well SPE plate and the 2 mg microelution 96-well SPE plate were both extremely effective at cleaning up our serum samples, the microelution format allowed us to produce ultra-concentrated samples without the need to dry down and reconstitute. By skipping this step we were able to save 30-45 minutes without losing any sensitivity. These time savings essentially doubled our productivity as the total cleanup time using the microelution format took approximately 30 minutes. In addition to the time savings, the microelution format is also amenable to small or limited sample volumes (as low as 10μ L) which provides an excellent solution in the event that our sample size decreases.



APPLICATIONS

Ordering Information

Strata®-X Microelution 96-Well SPE Plates				
Part No.	Description	Unit		
8M-S035-4GA	Strata-X-CW 33 µm Polymeric Weak Cation-Exchange Microelution 96-Well Plate, 2 mg/well	1/pk		
8M-S029-4GA	Strata-X-C 33 µm Polymeric Strong Cation-Exchange Microelution 96-Well Plate, 2 mg/well	1/pk		
8M-S100-4GA	Strata-X 33 µm Polymeric Reversed Phase Microelution 96-Well Plate, 2 mg/well	1/pk		
8M-S123-4GA	Strata-X-A 33 µm Polymeric Strong Anion-Exchange Microelution 96-Well Plate, 2 mg/well	1/pk		
8M-S038-4GA	Strata-X-AW 33 µm Polymeric Weak Anion-Exchange Microelution 96-Well Plate, 2 mg/well	1/pk		

Kinetex[®] Core-Shell HPLC/UHPLC Columns

5 µm Minibo	re Columns (mm)		SecurityGuard™ ULTRA Cartridges		dBore™ Columns (mm)		SecurityGuard ULTRA Cartridges [‡]
Phase	50 x 2.1	100 x 2.1	3/pk	Phase	50 x 3.0	100 x 3.0	3/pk
Biphenyl	00B-4627-AN	00D-4627-AN	AJ0-9209	Bipheny	00B-4627-Y0	00D-4627-Y0	AJ0-9208
			for 2.1 mm ID				for 3.0 mm ID
5 µm Analyti	cal Columns (mm)				SecurityGuard ULTRA Cartridges [‡]		
Phase	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk		
Biphenyl	00B-4627-E0	00D-4627-E0	00F-4627-E0	00G-4627-E0	AJ0-9207		
					for 4.6 mm ID		
2.6 µm Minib	oore Columns (mm)				SecurityGuard ULTRA Cartridges [‡]		
Phase	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk		
Biphenyl	00A-4622-AN	00B-4622-AN	00D-4622-AN	00F-4622-AN	AJ0-9209		
					for 2.1 mm ID		
2.6 µm MidB	ore Columns (mm)			SecurityGua ULTRA Cartrid			
Phase	50 x 3.0	100 x 3.0	150 x 3.0	3/pk			
Biphenyl	00B-4622-Y0	00D-4622-Y0	00F-4622-Y0	AJ0-9208			
				for 3.0 mm l	D		
2.6 µm Analytical Columns (mm)				SecurityGua ULTRA Cartrid			
Phase	50 x 4.6	100 x 4.6	150 x 4.6	3/pk			
Biphenyl	00B-4622-E0	00D-4622-E0	00F-4622-E0	AJ0-9207			
				for 4.6 mm l	D		
1.7 µm Minib	oore Columns (mm)			SecurityGua ULTRA Cartrid			
Phase	50 x 2.1	100 x 2.1	150 x 2.1	3/pk			
Biphenyl	00B-4628-AN	00D-4628-AN	00F-4628-AN	AJ0-9209			
		re holder. Part No.: AJ		for 2.1 mm II)		

[‡] SecurityGuard ULTRA Cartridges require holder, Part No.: AJ0-9000



If Phenomenex products in this technical note do not provide at least equivalent separations as compared to other products of the same phase, size and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

PI ICATIONS



Australia

t: +61 (0)2-9428-6444 f: +61 (0)2-9428-6445 auinfo@phenomenex.com

Austria

- t: +43 (0)1-319-1301 f: +43 (0)1-319-1300
- anfrage@phenomenex.com

Belaium

- t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) f: +31 (0)30-2383749 beinfo@phenomenex.com
- Canada
- t: +1 (800) 543-3681
- f: +1 (310) 328-7768 info@phenomenex.com

China

- t: +86 (0)20 2282-6668
- f: +86 (0)20 2809-8130 chinainfo@phenomenex.com

Denmark

- t: +45 4824 8048
- f: +45 4810 6265 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 f: +45 4810 6265 nordicinfo@phenomenex.com

- France t: +33 (0)1 30 09 21 10 f: +33 (0)1 30 09 21 11
- franceinfo@phenomenex.com

- **Germany** t: +49 (0)6021-58830-0
- f: +49 (0)6021-58830-11 anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400 f: +91 (0)40-3012 2411 indiainfo@phenomenex.com

Ireland

- t: +353 (0)1 247 5405
- f: +44 1625-501796 eireinfo@phenomenex.com

Italv

- t: +39 051 6327511 f: +39 051 6327555
- italiainfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Luxembourg t: +31 (0)30-2418700

f: +31 (0)30-2383749 nlinfo@phenomenex.com

Mexico t: 001-800-844-5226

- f: 001-310-328-7768
- tecnicomx@phenomenex.com

The Netherlands t: +31 (0)30-2418700 f: +31 (0)30-2383749

- nlinfo@phenomenex.com

New Zealand t: +64 (0)9-4780951

f: +64 (0)9-4780952 nzinfo@phenomenex.com

Norway t: +47 810 02 005

- f: +45 4810 6265
- nordicinfo@phenomenex.com

- Puerto Rico t: +1 (800) 541-HPLC f +1 (310) 328-7768
- info@phenomenex.com

Spain

- t: +34 91-413-8613 f: +34 91-413-2290
- espinfo@phenomenex.com

Sweden t: +46 (0)8 611 6950

- f: +45 4810 6265
- nordicinfo@phenomenex.com

United Kingdom

- t: +44 (0)1625-501367 f: +44 (0)1625-501796
- ukinfo@phenomenex.com

USA

- t: +1 (310) 212-0555
- f: +1 (310) 328-7768 info@phenomenex.com

All other countries Corporate Office USA t: +1 (310) 212-0555

- f: +1 (310) 328-7768
- info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex and Strata are registered trademarks and MidBore and SecurityGuard are trademarks of Phenomenex. API 5000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

© 2016 Phenomenex, Inc. All rights reserved.