


Etilefrine hydrochloride**Column:** N/A**Dimensions:** N/A**Order No:** N/A**Elution Type:** N/A**Eluent A:** Hexane**Eluent B:** 1,2-Dichloroethane**Eluent C:** Methanol**Eluent D:** Trifluoroacetic acid

Gradient	Step No.	Time (min)	Pct A	Pct B	Pct C	Pct D
Profile:	1	0	59.9	34.9	5	0.2

Flow Rate: 1 mL/min**Col. Temp.:** ambient**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 254 nm (ambient)**ANALYTES:**

- 1 Etilefrine hydrochloride
- 2 Etilefrine hydrochloride

