

Folic acid and metabolites using Luna 5u C18(2)

Column: Luna® 5 µm C18(2) 100 Å, LC Column 150 x 4.6 mm, Ea

Dimensions: 150 x 4.6 mm ID

Order No: 00F-4252-E0

Elution Type: Gradient

Gradient **Step No.** **Time (min)** **Pct A**

Profile:

Col. Temp.: ambient

Detection: SecurityGuard™ Guard Cartridges System x 200 µm I.D. (ambient)


- SecurityGuard Cartridges, C18 4 x 3.0mm, 10/Pk Part No.: AJ0-4287

- Holder Part No.: KJ0-4282

15736



Products used in this application:



ANALYTES:

- 1 unknown
- 2 unknown
- 3 5-Methyl tetrahydrofolic acid
- 4 Aminopterin
- 5 Folic acid
- 6 unknown
- 7 Methotrexate

