


Benzoylecgonine by MSD on a ZB-XLB

Column: N/A
Phase: N/A
Dimensions: N/A
Order No: N/A
Oven Profile: N/A

Carrier Gas: N/A
Injection: N/A
Detection: N/A

Analyst Note: Analytes are 24ppm derivatized with BSTFA
Sim ions:
(Quant ion*) D3-Benzoylecgonine (*243,259,364): Benzoylecgonine (*240,256,361)

16694

**ANALYTES:**

- 1** Benzoylecgonine-d3
- 2** Benzoylecgonine

