


Batch To Batch Reproducibility on bioZen Intact C4**Column:** N/A**Dimensions:** N/A**Order No:** N/A**Elution Type:** N/A**Eluent A:** 0.1% TFA in Water**Eluent B:** 0.1% TFA in ACN

Gradient	Step No.	Time (min)	Pct A	Pct B
Profile:	1	0	90	10
	2	0	70	30
	3	0	50	50
	4	0	90	10

Flow Rate: 1.5 mL/min**Col. Temp.:** 40 °C**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 280 nm (ambient)**ANALYTES:**

- 1** Ligh Chain
- 2** Heavy Chain
- 3** Intact mAb

