


Syringe Filter used in Separation of 23 PCB compounds on a 60m ZB-Dioxin GC column**Column:** Zebron[™] ZB-Dioxin, GC Cap. Column 60m, 0.25mm, 0.20µm, Ea**Phase:****Dimensions:** 60 meters x 0.25 mm x 0.2 µm**Order No:** 7KG-G045-10**Oven Profile:** °C for 2 min °C for 5 min**Carrier Gas:** ,**Injection:** :1 1 µL @ 280°C**Detection:** GC/MS (300°C)**Products used in this application:**

Syringe Filter used in Separation of 23 PCB compounds on a 60m ZB-Dioxin GC column

ANALYTES:

- 1 2,4',5-Trichlorobiphenyl
- 2 2,4,4'-Trichlorobiphenyl
- 3 2,3',4,6-Tetrachlorobiphenyl
- 4 2,2',5,5'-Tetrachlorobiphenyl
- 5 2,3',4',5-Tetrachlorobiphenyl
- 6 2,2',4,5,5'-Pentachlorobiphenyl
- 7 3,4,4',5-Tetrachlorobiphenyl
- 8 3,3',4,4'-Tetrachlorobiphenyl
- 9 2',3,4,4',5-Pentachlorobiphenyl
- 10 2,3',4,4',5-Pentachlorobiphenyl
- 11 2,2',4,4',5,5'-Hexachlorobiphenyl
- 12 2,3,3',4,4'-Pentachlorobiphenyl
- 13 2,3,3',4',5',6-Hexachlorobiphenyl
- 14 2,3,3',4',5,6-Hexachlorobiphenyl
- 15 2,2',3,4,4',5'-Hexachlorobiphenyl
- 16 3,3',4,4',5-Pentachlorobiphenyl
- 17 2,3',4,4',5,5'-Hexachlorobiphenyl
- 18 2,3,3',4,4',5-Hexachlorobiphenyl
- 19 2,3,3',4,4',5'-Hexachlorobiphenyl
- 20 2,2',3,4,4',5,5'-Heptachlorobiphenyl
- 21 3,3',4,4',5,5'-Hexachlorobiphenyl
- 22 2,2',3,3',4,4',5-Heptachlorobiphenyl
- 23 2,3,3',4,4',5,5'-Heptachlorobiphenyl

