


Gibberellin Homologues

Column: N/A**Dimensions:** N/A**Order No:** N/A**Elution Type:** N/A**Eluent A:** Acetic acid**Eluent B:** Water**Eluent C:** CH₃H

Gradient	Step No.	Time (min)	Pct A	Pct B	Pct C
Profile:	1	0	0.1	0.4	99.5

Flow Rate: 1.5 mL/min**Col. Temp.:** ambient**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 210 nm (ambient)**ANALYTES:****1** Gibberellin