


Hydrolyzed Dextran T-40

Column: N/A
Dimensions: N/A
Order No: N/A
Elution Type: N/A
Eluent A: Water
Gradient

Step No.	Time (min)	Pct A
1	0	100

Profile:
Flow Rate: 0.2 mL/min
Col. Temp.: 60 °C
Detection: Refractive Index (RI) @ nm (ambient)



ANALYTES:

- 1 Glucose
- 2 Isomaltose
- 3 Isomaltotriose
- 4 Isomaltotetraose
- 5 Isomaltopentaose
- 6 Isomaltohexaose
- 7 Isomaltoheptaose
- 8 Isomaltooctaose
- 9 Isomaltodecaose

