

Inorganic Anions

Column: Hamilton® 10 µm PRP-X100, LC Column 100 x 4.1 mm, Ea

Dimensions: 100 x 4.1 mm ID

Order No: CH0-0685

Elution Type: Isocratic

Eluent A: 4.0 mM p-Hydroxybenzoic acid (pH 8.5)

Eluent B: Methanol

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
1	0	97.5	2.5	

Flow Rate: 4 mL/min

Col. Temp.: ambient

Detection: UV-Vis Abs.-Variable Wave.(UV) @ 310 nm (ambient)

Analyst Note: SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, RP-1 4 x 3.0mm, 10/Pk Part No.: AJ0-5809

- Holder Part No.: KJ0-4282



Products used in this application:

Hamilton®

No Chromatogram Available

ANALYTES:

- 1 Fluoride
- 2 Carbonate
- 3 Chloride
- 4 Nitrite
- 5 Bromide
- 6 Nitrate
- 7 Phosphate
- 8 Sulfate

